

trūRapid™ HW/Lyme Test Product Bulletin

By Dr. Chinta M. Lamichhane B.V. Sc & AH (DVM), PhD, Frances Moore DVM DACVP, CQA (ASQ),
Dr. Haichen Song PhD, Xiuli Yang PhD, Maria Peluso-Lapsley, MBA

trūRapid Heartworm/Lyme test demonstrates excellent sensitivity and specificity for heartworm antigen and Lyme C6 antibody screening and detection.

Key Points

- Canine heartworm disease caused by the parasite *Dirofilaria immitis*, is a life-threatening disease.
- The prevalence of heartworm disease is increasing in North America and Europe.
- Lyme disease caused by the spirochete *Borrelia burgdorferi* (Bb) and related *Borrelia* strains results in lameness, fever, lethargy, anorexia and renal disease in dogs.
- The prevalence and geographic range of Lyme disease is increasing in North America, the United Kingdom (UK) and continental Europe.
- Heartworm antigen and Lyme/Bb C6 antibody assays in the trūRapid test are the recommended tests for screening and diagnosis of heartworm and Lyme disease.
- The trūRapid Heartworm/Lyme test has excellent sensitivity and specificity relative to reference methods and currently marketed commercial heartworm and Lyme tests.

Introduction

Dirofilaria immitis, the heartworm, is a parasitic roundworm that spreads from host to host through the bites of mosquitoes (Bowman D, et al, 2009). The parasite commonly resides in the pulmonary arterial system as well as the heart. A major effect on the health of infected animals is a manifestation of damage to the lung vessels and tissues (Dantas-Torres F, et al, 2013). Diagnosis of *Dirofilaria immitis* infection in dogs is based on many factors, including detection of an adult female worm antigen in serum, plasma, or anticoagulated whole blood. (<https://www.heartwormsociety.org/veterinary-resources/american-heartworm-society-guidelines>, 2024).

Lyme disease is a bacterial infection caused by *Borrelia burgdorferi* (Bb) sensu lato complex, bacterial spirochetes transmitted through the bite of an Ixodes sp. tick. In North America, *Borrelia burgdorferi* is the main cause of Lyme disease, while in Europe, *Borrelia afzelii* and *Borrelia garinii* are the primary causes. (Hengge, UR, et al 2003 ; Marques AR, et al 2021). Clinical signs of Lyme disease are most commonly reported in humans and in dogs. Asymptomatic infection is common in many domestic species. (Littman MP, et al, 2006; 2018; Springer A, et al, 2020). Wildlife, especially deer and mice, are reservoir hosts for Bb.

Heartworm

Heartworm Disease, caused by infection from *Dirofilaria immitis*, is a globally dispersed and clinically significant cardiopulmonary disease seen in canids and felids, as well as other genera (e.g. Mustelidae, Otariidae). Endemic areas in the United States have shown up to 45% prevalence, with tropical climates conferring the chance at even higher infection rates (Bowman D, et al, 2009) (Self SW, et al, 2019). Over the last 10 years, heartworm has spread from southern Europe into eastern and northeastern countries. (Morchon R et al, 2022).

Clinical signs of *D. immitis* infection in dogs include respiratory distress, epistaxis, hemoptysis, ascites, exercise intolerance, and anorexia. Despite the chance to observe those severe clinical signs, during early stages of infection and disease, dogs typically present asymptotically. Progression of disease tends to lead to worsening clinical signs and prognosis (<https://www.heartwormsociety.org/veterinary-resources/american-heartworm-society-guidelines>, 2024).

Medical and surgical management of canine heartworm disease can be difficult for a myriad reasons, including *D. immitis* life cycle (Table 1), side effects and efficacy of medical therapies, and clinical presentations in practice. Furthermore, even with complete resolution of the infection, some animals experience long term effects from the damage inflicted by the infection. Prevention and monitoring remain key aspects of eradication programs, regardless of levels of endemicity (<https://www.heartwormsociety.org/veterinary-resources/american-heartworm-society-guidelines>, 2024).

Table 1. *Dirofilaria immitis* Life Cycle

<i>D. Immitis</i> Life Stage	L1, L2, L3	L4	S5 (immature adult)	Adult
Stage Complete	10–14 Days	45–65 Days	4–5 Months	6– months (total)
Location	Mosquito	Host tissue	Host bloodstream	Host bloodstream

Point-of-care heartworm antigen tests can detect *D. immitis* in whole blood, serum or plasma, and give veterinary professionals vital tools in diagnosing and monitoring canine heartworm infection and disease. Antigen tests are developed to recognize

antigen of adult female heartworms, which can complicate diagnosis in the presence of low worm burdens, infections with only larvae or immature stages of the worms, or all-male infections. Sensitivity and specificity remain high for antigen tests broadly, with as much as 99% sensitivity and 100% specificity, depending on test type and worm burden (<https://www.heartwormsociety.org/veterinary-resources/american-heartworm-society-guidelines>, 2024)

Lyme/Borrelia sp.

The prevalence and geographical range of Lyme disease is increasing although precise estimates are difficult due to variation in case definition and surveillance methods. (Rizzoli A, et al, 2011; Dantes – Torres F et al, 2012; Bajer A, et al, 2014; Michelet L, et al, 2014; Bjurman NK, et al, 2016; Cairns V, et al, 2019; Tulloch JSP et al, 2019; Springer A et al, 2020; Marques, AR, et al, 2021; Miro G et al, 2022; Kugeler KJ et al, 2022; Bajer A et al, 2022; <https://capcvet.org/>; <https://www.canada.ca/en/public-health/services/diseases/lyme-disease/surveillance-lyme-disease.html#a1>) The Mid-Atlantic, Northeast, Upper Midwest regions of the United States have the greatest prevalence of Lyme disease in the United States. (Kugeler KJ, et al, 2022) In Canada, most cases are reported in Ontario, Quebec and Nova Scotia. (<https://www.canada.ca/en/public-health/services/publications/diseases-conditions/lyme-disease-surveillance-canada-annual-edition-2022.html>) In the UK, the majority of cases are reported from Scotland, South and Southwest England. (Cairns V, et al, 2019; Tulloch JSP et al, 2019) In continental Europe, the prevalence of Lyme Disease is greatest in Scandinavia, Baltic States, Austria, Czech Republic, Germany and Slovenia. (Marques AR, et al, 2021) Expansion of the range of the Ixodes sp. tick vectors associated with climate change, encroachment of suburban into rural areas and increased human activities in tick environments are important factors in the increased number of Bb infections. (Littman MP, et al, 2018)

Diagnosis of Lyme disease depends on demonstration of Bb antibodies in the blood of dogs with characteristic clinical signs. The clinical signs and syndromes most associated with Lyme disease are lameness, fever, lethargy, anorexia, lymph node enlargement and glomerulonephritis. (Littman MP, 2006, 2018) However, the majority of Bb infections in animals are asymptomatic which presents challenges in confirming a diagnosis. Serology for the Bb C6 antibody is the recommended assay for supporting a diagnosis of Lyme disease and eliminates interference resulting from vaccine associated antibodies. (Littman MP, et al, 2018) The C6 peptide is highly specific and antigenically conserved among various species of *B. burgdorferi* sensu lato complex. The C6 peptide derived from *B. burgdorferi* detects the infection regardless of infecting strains. (Liang FT, et al., 2000. Tjernberg I, et al 2007; Christova I et al, 2013; Rouhiainen M et al, 2021). Levels of C6 antibody decrease in infected animals during treatment. The trūRapid Heartworm/Lyme test detects the C6 antibody. Routine serologic screening of dogs is recommended in areas endemic for Lyme disease, followed by screening seropositive dogs for proteinuria. (Littman MP, 2018) The presence of Bb antibody indicates tick exposure and a risk for not only Lyme disease but also for a variety of other tick-borne diseases. (Dantes – Torres F, et al, 2012; Michelet L, et al, 2014; Bajer A, et al, 2022; Miro G et al, 2022).


Confirming a diagnosis of Lyme disease in a symptomatic dog requires demonstration of C6 antibody and rule out of other tick-borne diseases or other potential causes of the presenting clinical signs. Initiation of measures to limit tick exposure is warranted in any dog seropositive for Bb.

trūRapid HW/Lyme technology

trūRapid HW/Lyme point of care (POC) test strips are sandwich-immunoassays intended for the visual qualitative detection of both heartworm antigen (*Dirofilaria immitis*) and canine *Borrelia burgdorferi* antibody in canine serum, plasma or anticoagulated whole blood on a single nitrocellulose membrane-based test strip.

The assays are set up as follows: two antibodies are used in the test strip for the detection of the female *D. immitis* antigen. One antibody is coupled on gold nanoparticles and a second, high specificity antibody is immobilized in the heartworm test line area (H). A combination of a protein and a highly purified Lyme C6 peptide are used in the test strip for *Borrelia burgdorferi* antibody detection. The purified peptide is coupled on gold nanoparticles and the protein is immobilized in the Lyme test line area (L). The control gold nanoparticles bind to the antigen, a protein immobilized in the control line area, forming a control line.

When the sample, followed by the sample buffer, is put into the



sample wells, it will be absorbed into the absorbing pad of the test strip. The fluid mixes with the gold labelled antibodies and peptides in the conjugate pad. Due to the capillary action of the test strip, the fluid will run up over the test strip, crossing the two test line regions (H and L), then subsequently the control line region. The control line should always appear, signifying correct functioning of the test

If the patient sample contains the antigen of even a single adult female heartworm, a line will appear in the heartworm test line region (H). This happens by building a sandwich between the gold labelled antibodies from the conjugate pad, the antigen from the sample, and the immobilized antibody of the heartworm test line region. Even a faint line is considered positive for heartworm infection.

If no antigen is in the sample, the gold labelled antibodies cannot connect to the immobilized antibodies in the heartworm test line region, and therefore no heartworm test line appears.

If the patient sample contains the *Borrelia burgdorferi* C6 antibody, a line will appear in the Lyme test line region (L). This happens by building a sandwich between the gold labelled peptide from the conjugate pad, the Lyme antibody from the sample, and the immobilized protein of the Lyme test line region. If the test is positive for *Borrelia burgdorferi* C6 antibodies, it means that the animal has antibodies from natural exposure to *Borrelia burgdorferi* (Littman MP, et al, 2018).

If no *Borrelia burgdorferi* C6 antibody is in the sample, the gold labelled peptide cannot connect to the immobilized protein in the Lyme test line region and therefore no Lyme test line appears.

Material and Methods

trūRapid HW/Lyme test strips were tested in this performance evaluation. A total of 270 and 128 samples from naturally infected dogs were tested in two comparison studies comparing the trūRapid HW/Lyme results to other methods. The heartworm performance data was generated by testing samples using a commercial enzyme linked immunosorbent assay (ELISA) and a rapid antigen test. The comparisons for Lyme were generated by testing against an indirect fluorescent antibody assay (IFA) and a rapid antibody test.

Results

Sensitivity and Specificity

Study 1: The sensitivity and specificity of the test kit were evaluated by testing a total of 270 dog serum samples with trūRapid HW/Lyme tests. Sensitivity and specificity were determined comparing DiroCHEK® HW Antigen test, as shown with a 2x2 analysis table (Table 1 and 2). All data is from naturally infected animals.

Table 1: Comparison of heartworm test results - trūRapid versus DiroCheck

trūRapid Heartworm Test	DiroChek® Heartworm Antigen Test		Total
	Detected	Undetected	
Detected	150	2	152
Undetected	0	118	118
Total	150	120	270
Sensitivity	100% (95% CI: 97.6%–100%)		
Specificity	98.3% (95% CI: 94.1%–99.8%)		

Table 2: Comparison of Lyme test results—trūRapid versus IFA

trūRapid Lyme Test	Lyme IFA ¹		Total
	Detected	Undetected	
Detected	103	3	106
Undetected	1	163	164
Total	104	166	270
Sensitivity	99.0% (95% CI: 94.8%–100%)		
Specificity	98.2% (95% CI: 94.8%–99.6%)		

¹VMRD - SLD-IFA-LD test.

Study 2: The performance of the trūRapid HW/Lyme test kit was evaluated by testing a total of 128 dog serum samples for both heartworm and Lyme. For heartworm, an additional 150 positive samples were tested. The sensitivity and specificity were compared with IDEXX SNAP 4DX®, with a 2x2 analysis table. All data is from naturally infected animals.

Table 3: Comparison of Heartworm test results - trūRapid HW/Lyme test versus SNAP 4Dx

trūRapid HW/Lyme Test (Heartworm results)	IDEXX SNAP 4DX (Heartworm Results)		
	Detected	Undetected	Total
Detected	149	1	150
Undetected	1	127	128
Total	150	128	278
Sensitivity	99.3% (95% CI: 96.3%–100%)		
Specificity	99.2% (95% CI: 95.7%–100%)		

Table 4: Comparison of Lyme test results - trūRapid HW/Lyme test versus SNAP 4Dx

trūRapid HW/Lyme Test (Lyme results)	IDEXX SNAP 4DX (Lyme Results)		
	Detected	Undetected	Total
Detected	49	2	51
Undetected	1	76	77
Total	50	78	128
Sensitivity	98.0% (95% CI: 89.4%–100%)		
Specificity	97.4% (95% CI: 91.0%–99.7%)		

Conclusion

The study demonstrates Antech trūRapid HW/Lyme test delivers excellent sensitivity and specificity for heartworm antigen and Lyme C6 antibody when compared to reference methods and in-house, rapid test methods. It is reliable test that provides rapid results for screening and detection.

Acknowledgements

The authors would like to thank Dr. Yanjin Zhang, Associate Professor, Department of Veterinary Medicine, VA-MD Regional College of Veterinary Medicine College Park Maryland USA and Dr. Jaswinder Sain DVM, Cascades Pet Hospital, VA for their time and dedication to this process. They provided invaluable assistance to study design, implementation, analysis, and publication.

References (Heartworm)

1. American Heartworm Society Canine Guidelines for the Prevention, Diagnosis and Management of Heartworm (*Dirofilaria immitis*) Infection in Dogs. Revised 2024 (<https://www.heartwormsociety.org/veterinary-resources/american-heartworm-society-guidelines>)
2. Bowman D, Little SE, Lorentzen L, Shields J, Sullivan MP, Carlin EP. Prevalence and geographic distribution of *Dirofilaria immitis*, *Borrelia burgdorferi*, Ehrlichia canis and Anaplasma phagocytophilium in dogs in the United States: result of a national clinic-based serologic survey. *Vet Parasitol.* 2009; 9; 160:138-48.
3. Dantas-Torres F, Otranto D. *Dirofilaria immitis*? Parasit Vectors. 2013 Oct 2;6(1):288.
4. Morchón R, Montoya-Alonso JA, Rodríguez-Escolar I, Carretón E. What Has Happened to Heartworm Disease in Europe in the Last 10 Years? *Pathogens.* 2022 Sep 13;11(9):1042. doi: 10.3390/pathogens11091042. PMID: 36145474; PMCID: PMC9503846.
5. Rohrbach BW, Patton S. Effects of diagnostic test accuracy and treatment efficacy on the occurrence of suspected failure of heartworm prophylaxis in dogs. *J Vet Intern Med.* 2013; 27(4):791-7.
6. Self, S.W., Pulaski, C.N., McMahan, C.S. et al. Regional and local temporal trends in the prevalence of canine heartworm infection in the contiguous United States: 2012–2018. *Parasites Vectors* 12, 380 (2019). <https://doi.org/10.1186/s13071-019-3633-2>

References (Lyme)

1. Bajer A, Kowalec M, Levytska VA, Mierzejewska EJ, Alsarraf M, Poliukhovych V, Rodo A, Wężyk D, Dwużnik-Szarek D. Tick-Borne Pathogens, Babesia spp. and *Borrelia burgdorferi* s.l., in Sled and Companion Dogs from Central and North-Eastern Europe. *Pathogens.* 2022 Apr 21;11(5):499. doi: 10.3390/pathogens11050499. PMID: 35631020; PMCID: PMC9144822.
2. Bajer A, Mierzejewska EJ, Rodo A, Bednarska M, Kowalec M, Welc-Falęciak R.
3. The risk of vector-borne infections in sled dogs associated with existing and new endemic areas in Poland: Part 1: A population study on sled dogs during the racing season. *Vet Parasitol.* 2014 May 28;202(3-4):276-86. doi: 10.1016/j.vetpar.2013.12.033. Epub 2014 Jan 20. PMID: 24491396.
4. Barbour AG, Zückert WR. Genome sequencing. New tricks of tick-borne pathogen. *Nature.* 1997 Dec 11;390(6660):553, 555. doi: 10.1038/37475. PMID: 9403678.
5. Bjurman NK, Bradet G, Lloyd VK. Lyme disease risk in dogs in New Brunswick. *Can Vet J.* 2016 Sep;57(9):981-4. PMID: 27587892; PMCID: PMC4982571. *BMJ Open.* 2019 Jul 31;9(7):e028064. doi: 10.1136/bmjopen-2018-028064. PMID: 31362976; PMCID: PMC6677960.
6. Cairns V, Wallenhorst C, Rietbrock S, Martinez C. Incidence of Lyme disease in the UK: a population-based cohort study. *BMJ Open.* 2019 Jul 31;9(7):e025916. doi: 10.1136/bmjopen-2018-025916. PMID: 31362975; PMCID: PMC6677952.
7. Christova, Iva & Trifonova, Iva & Gladnishka, Teodora & Taseva, Evgenia & Ivanova, Veda & Rusimova, Dobrinka. (2014). C6 Peptides from *Borrelia burgdorferi* Sensu Stricto, *Borrelia Afzelii* and *Borrelia Garinii* as Antigens for Serological Diagnosis of Lyme Borreliosis. *Biotechnology & Biotechnological Equipment.* 27. 3540-3542. 10.5504/BBEQ.2012.0100.
8. Dantas-Torres F, Chomel BB, Otranto D. Ticks and tick-borne diseases: a One Health perspective. *Trends Parasitol.* 2012 Oct;28(10):437-46. doi: 10.1016/j.pt.2012.07.003. Epub 2012 Aug 17. Erratum in: *Trends Parasitol.* 2013 Oct;29(10):516. PMID: 22902521.
9. Hengge UR, Tannapfel A, Tying SK, Erbel R, Arendt G, Ruzicka T. Lyme borreliosis. *Lancet Infect Dis.* 2003 Aug;3(8):489-500. doi: 10.1016/s1473-3099(03)00722-9. Erratum in: *Lancet Infect Dis.* 2003 Dec;3(12):815. PMID: 12901891.

11. <https://capcvet.org/>
12. <https://www.canada.ca/en/public-health/services/diseases/lyme-disease/surveillance-lyme-disease.html#a1>
13. Kugeler KJ, Earley A, Mead PS, Hinckley AF. Surveillance for Lyme Disease After Implementation of a Revised Case Definition — United States, 2022. *MMWR Morb Mortal Wkly Rep* 2024;73:118–123. DOI: <http://dx.doi.org/10.15585/mmwr.mm7306a1>
14. Liang FT, Aberer E, Cinco M, Gern L, Hu CM, Lobet YN, Ruscio M, Voet PE Jr, Weynants VE, Philipp MT. Antigenic conservation of an immunodominant invariable region of the VlsE lipoprotein among European pathogenic genospecies of *Borrelia burgdorferi* SL. *J Infect Dis.* 2000 Nov;182(5):1455-62. doi: 10.1086/315862. Epub 2000 Oct 9. PMID: 11023468.
15. Littman MP, Gerber B, Goldstein RE, Labato MA, Lappin MR, Moore GE. ACVIM consensus update on Lyme borreliosis in dogs and cats. *J Vet Intern Med.* 2018 May;32(3):887-903. doi: 10.1111/jvim.15085. Epub 2018 Mar 22. PMID: 29566442; PMCID: PMC5980284.
16. Littman MP, Goldstein RE, Labato MA, Lappin MR, Moore GE. ACVIM small animal consensus statement on Lyme disease in dogs: diagnosis, treatment, and prevention. *J Vet Intern Med.* 2006 Mar-Apr;20(2):422-34. doi: 10.1892/0891-6640(2006)20[422:asacso]2.0.co;2. PMID: 16594606.
17. Marques AR, Strle F, Wormser GP. Comparison of Lyme Disease in the United States and Europe. *Emerg Infect Dis.* 2021 Aug;27(8):2017-2024. doi: 10.3201/eid2708.204763. PMID: 34286689; PMCID: PMC8314816.
18. Michelet L, Delannoy S, Devillers E, Umhang G, Aspan A, Juremalm M, Chirico J, van der Wal FJ, Sprong H, Boye Pihl TP, Klitgaard K, Bødker R, Fach P, Moutailler S. High-throughput screening of tick-borne pathogens in Europe. *Front Cell Infect Microbiol.* 2014 Jul 29;4:103. doi: 10.3389/fcimb.2014.00103. PMID: 25120960; PMCID: PMC4114295.
19. Miró G, Wright I, Michael H, Burton W, Hegarty E, Rodón J, Buch J, Pantchev N, von Samson-Himmelstjerna G. Seropositivity of main vector-borne pathogens in dogs across Europe. *Parasit Vectors.* 2022 Jun 6;15(1):189. doi: 10.1186/s13071-022-05316-5. PMID: 35668469; PMCID: PMC9169295.
20. Rouhiainen M, Pietikäinen A, Kortela E, Kanerva MJ, Oksi J, Hytönen J. C6 peptide enzyme immunoassay in Lyme borreliosis serology. *J Microbiol Methods.* 2021 Jan;180:106122. doi: 10.1016/j.mimet.2020.106122. Epub 2020 Dec 14. PMID: 33326821.
21. Speck S, Reiner B, Streich WJ, Reusch C, Wittenbrink MM. Canine borreliosis: a laboratory diagnostic trial. *Vet Microbiol.* 2007 Feb 25;120(1-2):132-41. doi: 10.1016/j.vetmic.2006.10.017. Epub 2006 Oct 21. PMID: 17101241.
22. Springer A, Glass A, Topp AK, Strube C. Zoonotic Tick-Borne Pathogens in Temperate and Cold Regions of Europe—A Review on the Prevalence in Domestic Animals. *Front Vet Sci.* 2020 Dec 10;7:604910. doi: 10.3389/fvets.2020.604910. PMID: 33363242; PMCID: PMC7758354.
23. Tjernberg, I., Krüger, G. & Eliasson, I. C6 peptide ELISA test in the serodiagnosis of Lyme borreliosis in Sweden. *Eur J Clin Microbiol Infect Dis* 26, 37–42 (2007). <https://doi.org/10.1007/s10096-006-0239-3>
24. Tulloch JSP, Semper AE, Brooks TJG, Russell K, Halsby KD, Christley RM, Radford AD, Vivancos R, Warner JC. The demographics and geographic distribution of laboratory-confirmed Lyme disease cases in England and Wales (2013–2016): an ecological study. *BMJ Open.* 2019 Jul 31;9(7):e028064. doi: 10.1136/bmjopen-2018-028064. PMID: 31362976; PMCID: PMC6677960.