

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT ME0125A

 2025-03-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0125A-1	LOT	ME0125A-2	LOT	ME0125A-3
element HT5+	Ret# x10 ³ /μL	76.4	± 64.4	214.3	± 99.0	189.9	± 99.0
QC Mode	Ret%	1.29	± 1.12	4.69	± 2.00	8.15	± 3.50
	IRF %	12.0	± 15.0	12.9	± 13.0	9.6	± 25.0
	LFR %	88.0	± 15.0	87.1	± 15.0	90.4	± 25.0
	MFR %	10.7	± 12.0	11.5	± 10.0	7.0	± 22.0
	HFR %	1.3	± 10.0	1.4	± 10.0	2.6	± 10.0
	RHE pg	24.6	± 12.0	22.5	± 12.0	19.6	± 12.0
	IPF %	7.5	± 15.0	7.6	± 15.0	7.6	± 15.0

Before using, refer to the instruction sheet for mixing directions.