

# BC-RET

## HEMATOLOGY CONTROLS

**CONTROL**

ASSAY VALUES AND EXPECTED RANGES

**LOT** ME1124A



2025-01-10

Instrument	Parameter	Low		Normal		High	
		<b>LOT</b>	ME1124A-1	<b>LOT</b>	ME1124A-2	<b>LOT</b>	ME1124A-3
element HT5+	Ret# x10 <sup>3</sup> /μL	77.5	± 64.4	226.6	± 99.0	197.0	± 99.0
QC Mode	Ret%	1.35	± 1.12	4.81	± 2.00	8.21	± 3.50
	IRF %	13.5	± 15.0	10.9	± 13.0	9.0	± 25.0
	LFR %	86.5	± 15.0	89.1	± 15.0	91.0	± 25.0
	MFR %	12.4	± 12.0	9.0	± 10.0	6.5	± 22.0
	HFR %	1.1	± 10.0	1.9	± 10.0	2.5	± 10.0
	RHE pg	26.0	± 12.0	23.5	± 12.0	20.0	± 12.0
	IPF %	7.4	± 15.0	7.5	± 15.0	7.5	± 15.0

Before using, refer to the instruction sheet for mixing directions.